## Cabinet – Type 2 – 25-604535-4/-5



## System 25

All components are RoHS-compatible.

## Design:

- Work cabinet for placing on a table/plane surface
- Dimensions of the base: 635 x 480 mm
- · Total height: 350 mm
- Made of vacuum moulded transparent PETG
- Provided with a pane made of safety glass
- Pane corners made of PP in the colours white (-5) or red (-4)
- Recommended airflow: 117 m<sup>3</sup>/h
- Optional flange for extraction available



## **Working area**

