



System 25

All components are RoHS-compatible.

Design:

- Work cabinet for placing on a table/plane surface
- Dimensions of the base: 1000 x 600 mm
- Height of side panels: 200 mm (total height 410 mm), art. no. 25-106020-4/-5
500 mm (total height 710 mm), art. no. 25-106050-4/-5
750 mm (total height 960 mm), art. no. 25-106075-4/-5
- Vacuum moulded top made of transparent PETG with a pane of tempered glass
- Pane corners made of PP in the colours white (-5) or red (-4)
- Sides made of transparent PETG
- Corners made of anodized aluminium
- Bottom edge provided with a non-skid rubber edge
- Recommended airflow: 220, 550 or 850 m³/h dependent on the height of the side panels
- Optional flange for extraction available
- NB! Can be delivered with doors



Working area

