



System 25

All components are RoHS-compatible.

Design:

- Work cabinet for placing on a table/plane surface
- Dimensions of the base: 420 x 320 mm
- Height of side panels: 400 (total height: 515 mm)
- Vacuum moulded top made of transparent PETG inclusive of 2 stub pipes
- Sides made of steel with a white polyester powder-coating on the surface
- Flanges, stub pipes and joint with damper Ø75 mm
- Recommended airflow: 204 m³/h



NB! Reductions to be ordered separately

Working area

