



## System 50 Flex - Aluminium (Flex AL)

All components RoHS-compatible.

### Design:

- Extraction arm for mounting on table
- Working range 900 mm
- Recommended airflow 45 to 85 m<sup>3</sup>/h
- Diameter of tube: Ø50 mm
- Tubes and extractor tube made of anodized Aluminium - surface 10 µ - in combination with a flexible hose
- Increased options of positioning due to the flexible hose
- Joints made of shatterproof, chemical resistant polypropylene (PP)
- With funnel of the extractor tube in order to increase the capture efficiency
- All O-rings maintenance free
- Integrated valve. When in an open position out of the flow path to maintain the maximum amount of flow
- All threaded stays, springs and thumbscrews made of acid-proof stainless steel (AISI 316L)
- Acid-proof POP<sup>®</sup> rivets for increased durability in aggressive environments
- Colour of joints/flange: white (-5), red (-4) or black (-050)
- Dismantling of the arm without tools for cleaning



Links: [Pressure drop chart](#), [Mounting](#), [User's manual](#), [Suction pen](#)

## Working area

Units in mm

