



## System 100 - Aluminium (AL)

- Round hood recommended for light fumes, gases and over small open vessels
- Diameter of hood: Ø385 mm
- Increased stability when moving the hood due to reinforced rim of the hood
- Increased efficiency at an angled position
- Hood made of transparent PETG, resistant to solvents
- Connection tube made of anodized aluminium - surface 10 µ
- The transparent hood assists in keeping a good view of the work place



## Capture efficiency

For this product there are no measurements.

In order to achieve optimal position, please compare with no. [1-10050](#).

## Drawing

